



Distillation apparatus set

Αρ. προϊόντος:

2265300

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής

Fast set-up for sample pretreatment.

Fast set-up for sample pretreatment. This Distillation set includes EPA approved glassware ONLY. Easy to set up because you only receive the parts you need. Illustrated, step-by-step procedures simplify distillation.

*Heater and support not included.

Distillation set includes all the glassware you need to satisfy EPA distillation requirements

Easy to set up - you only receive the parts you need

Illustrated, step-by-step procedures simplify distillation

Προδιαγραφές

Όγκος/Μέγεθος πακέτου :

Each 1

Τι περιέχει το καλάθι:

Distillation flask, J-tube, condenser, drip tube, erlenmeyer flask, and non-Hg thermometer

Τι περιέχει το καλάθι

Distillation flask, J-tube, condenser, drip tube, erlenmeyer flask, and non-Hg thermometer