



REF921 Reference electrode, calomel w/ LiCl solution, screw cap (Radiometer Analytical)

Αρ. προϊόντος:

E21M007

ΚΑΤΑΡΓΗΜΕΝΟ ΣΤΟΙΧΕΙΟ

Αυτό το στοιχείο δεν είναι πλέον διαθέσιμο.

Reference electrode for non-aqueous media.

The right probe for your application.

To meet your analysis needs, Radiometer Analytical reference electrodes are available with a variety of reference systems and liquid junctions. There is a Radiometer electrode for every application.