



Be Right™



Sension+ 9660C Combination Calcium ISE

Αρ. προϊόντος:

LZW9660C.97.002

EUR Τιμή:

Επικοινωνήστε μαζί μας

Διαθέσιμο

Utilizes a unique solid-state sensor technology that eliminates costly membrane replacements

Traditional PVC-membrane ISEs require membrane replacement every 2-3 months. The 9660C ISE's unique solid-state sensor technology eliminates the need for frequent membrane replacement by using a solid gel ion exchange behind the ion-selective membrane versus the liquid ion exchange used in most ISEs. This design allows users to get up to 2 years of life from their ISE without membrane replacement – saving the time and costly expense of frequently exchanging membranes

Requires virtually no maintenance

The solid-state sensor design allows for DRY storage of the ISE without a shelf life or membrane replacement.

Provides fast, stable, and accurate response in a variety of sample types

Προδιαγραφές

Accuracy:	±0.02 mV or 0.05%, whichever is greater (application dependent)
Electrode junction :	Porous PTFE Annular Ring
Reference system:	Ag/AgCl (double junction)
Απαιτείτε ISA:	3M KCl
Διάλυμα πλήρωσης :	Non-refillable driTEK gel
Διαστάσεις (BxM):	12 mm x 120 mm
Εγγύηση :	6 μήνες
Ειδικά χαρακτηριστικά :	Combination electrode. Maintenance free: No electrolyte or spare membranes needed
Εύρος θερμοκρασίας :	0 - 50 °C
Εύρος μέτρησης :	0.02 - 4000 mg/L Ca ²⁺
Θερμική αντίσταση :	-
Μήκος :	120 mm
Μήκος καλωδίου :	1 m
Παράμετρος :	Calcium
Σύνδεσμος :	BNC Connector
Τύπος :	Εργαστηριακό ISE συνδυαστικό
Τύπος αισθητηρίου :	Solid-state PVC membrane
Τύπος ηλεκτροδίων :	Laboratory Combination ISE

Απαιτούμενα παρελκόμενα

- Διάλυμα αποθήκευσης ηλεκτροδίου pH, 500 mL (Item 2756549)