



Beaker, polypropylene, low form, 50 mL

Αρ. προϊόντος:

108041

EUR Τιμή:

Επικοινωνήστε μαζί μας

Διαθέσιμο

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes ≥ 250 mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Προδιαγραφές

Διάμετρος, εξωτερική :	51.1 mm with spout
Μονάδες μέτρησης:	1 each
Όγκος/Μέγεθος πακέτου :	1 each
Σώμα υλικού:	Polypropylene
Υλικό :	Polypropylene, Low Form
Χωρητικότητα :	50 mL