



Buret, Acrylic, 50 mL

Αρ. προϊόντος:

1405941

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής

Durable, crystal-clear acrylic body is individually calibrated to meet the accuracy requirements of ASTM E287. Leak-proof stopcock with self-lubricating PTFE TFE plug. Excellent for general lab or field use in schools, industrial, municipal or food applications. Unaffected by dilute mineral acids and bases (except hydrofluoric acid and ammonium hydroxide). Can contain liquids within a pH range of 1 - 14 and up to 5% wt/wt of acids or 30% bases can be used. Not designed for exposure to alcohol or organic solvents.

Προδιαγραφές

Graduations:	0.1 mL
Ειδικά χαρακτηριστικά :	Buret
Σώμα υλικού:	Acrylic
Χωρητικότητα :	50 mL