



Manganese test kit, model MN-5

Αρ. προϊόντος:

146700

Επικίνδυνο

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής



Manganese causes a bitter taste in water, and at concentrations above 0.1 mg/L, it causes objectionable stains on laundry and plumbing fixtures. Manganese control also is important in the beverage, paper, and textile industries, and in dye production and food processing. An economical solution, this kit uses the simple, accurate color disc method.

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Προδιαγραφές

Parameter:	Manganese
Smallest increment:	0.1
Αριθμός δοκιμών :	100
Εύρος μέτρησης :	0 - 3.0 mg/L Mn
Κατηγορία:	Χρωματικός δίσκος
Μοντέλο :	MN-5
Μοντέλο σετ μέτρησης :	MN-5
Όνομα μεθόδου :	Color disc/Cold periodate
Παράμετρος :	Manganese - As Mn
Τυπικό εύρος μέτρησης:	0 - 3.0 mg/L Mn
Τύπος θήκης :	D