



Laboratory Bottles, Polypropylene, Autoclavable, Wide Mouth, 1 L, 6/pk

Αρ. προϊόντος: 1472463

EUR Τιμή: Επικοινωνήστε μαζί μας
Δεν εμφανίζεται η ημερομηνία αποστολής

Polypropylene bottles suitable for sample collection to chemical storage. Matched polypropylene caps have a valve-seal design for a perfect, leakproof fit. A protruding ring under the cap holds shrink-wrapping securely in place. Autoclavable 121 °C (Do not autoclave with caps in place). Available in Narrow Mouth or Wide Mouth.

Προδιαγραφές

Body Material:	Plastic
Package size UOM:	6pk
Όγκος/Μέγεθος πακέτου :	6 /pk
Υλικό :	PP
Χωρητικότητα :	1000 mL