



## Chromium color disc test kit, model CH-8

Αρ. προϊόντος:

183400

Επικίνδυνο

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής



Widely used as a corrosion inhibitor in cooling towers and engine cooling systems, chromium is also used in surface finishing, leather tanning, and other industries. The most widely used form, hexavalent chromium, is very toxic and strict limits are placed on discharges. Utilizing the color disc method to determine chromium (Range: 0.1-1.5 mg/L), this kit contains 100 tests.

**Simple, low-cost color disc method is the most accurate visual test kit method**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

---

### Προδιαγραφές

Parameter:	Chromium Hexavalent
Smallest increment:	0.1
Αριθμός δοκιμών :	100
Εύρος μέτρησης :	0 - 1,3 Cr <sup>6+</sup>
Κατηγορία:	Χρωματικός δίσκος
Μοντέλο :	CH-8
Μοντέλο σετ μέτρησης :	CH-8
Όνομα μεθόδου :	Color disc/ Diphenylcarbazide
Παράμετρος :	Chromium, low range - As hexavalent Cr
Τύπος θήκης :	D