



Hydrazine color disc test kit, model HY-2

Αρ. προϊόντος:

184900

Επικίνδυνο

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής



Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. Utilizing the color disc method with a range of 0.02-1 mg/L, this kit contains 300 tests.

Measures Hydrazine, a common oxygen scavenger in boiler water, using the simple color disc method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Προδιαγραφές

Parameter:	Detergents
Smallest increment:	0
Αριθμός δοκιμών :	300
Εύρος μέτρησης :	0 - 1 mg/L N_2H_4
Κατηγορία:	Χρωματικός δίσκος
Μοντέλο :	HY-2
Μοντέλο σετ μέτρησης :	HY-2
Όνομα μεθόδου :	Color disc/
	p-Dimethylamino-
	benzaldehyde
Παράμετρος :	Hydrazine - As N_2H_4
Τυπικό εύρος μέτρησης:	0 - 1 mg/L N_2H_4
Τύπος θήκης :	D