



## Test kit tannin/lignin

**Αρ. προϊόντος:**

**193701**

**Επικίνδυνο**

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής



A product of decomposed plant material, tannin occurs in natural waters and is used for boiler and cooling tower water treatment, dyeing, ink manufacture, tanning, and paper sizing. Lignin, also a plant product, is found in paper pulp manufacturing waste. Tannin and lignin are analyzed together and reported as mg/L tannic acid. This kit utilizes the simple, low-cost color disc method - the most accurate visual test kit method.

**Tannin and lignin are analyzed together and reported as mg/L tannic acid.**

**Simple, low-cost color disc method is the most accurate visual test kit method**

**Achieve fast results with an easy-to-carry test kit**

## Προδιαγραφές

Parameter:	Tannin/Lignin
Smallest increment:	0.5
Αριθμός δοκιμών :	5
Εύρος μέτρησης :	100
Εύρος μέτρησης (2):	0 - 15 mg/L as Tannic Acid
Κατηγορία:	0 - 150 mg/L
Μοντέλο :	Χρωματικός δίσκος
Μοντέλο σετ μέτρησης :	Tannin/Lignin
Όνομα μεθόδου :	TA-3
Παράμετρος :	Color disc/Tyrosine
Τυπικό εύρος μέτρησης:	Tannin/Lignin - As tannic acid
Τύπος θήκης :	0 - 15 mg/L as Tannic Acid; 0 - 150 mg/L
	R