



Be Right™



## Test kit, hardness & iron color disc, model HA-77

Αρ. προϊόντος:

202300

Επικίνδυνο

EUR Τιμή:

Επικοινωνήστε μαζί μας

Τηλεφωνήστε για την ημερομηνία αποστολής



Hardness & Iron Water Quality Test Kit is the kit to use if you need better accuracy than that provided by a test cube. Utilizing drop count titration and color disc colorimetry, this kit contains 100 tests.

**The kit to use for more resolution in a visual test**

**Uses simple, economical drop count titration and color disc methods**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

---

### Προδιαγραφές

Εφαρμογή:

Water Quality

Κατηγορία:

Χρωματικός δίσκος