



Test kit, acidity, model AC-DT

Αρ. προϊόντος:

2064000

Επικίνδυνο

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής



Acidity is usually caused by the presence of mineral acids, salts of strong acids, or free carbon dioxide. In surface or ground waters, acidity may form natural substances or industrial pollution. Excessive levels cause corrosion and can be detrimental to fish. Utilizing the Hach Digital Titrator, the Model AC-DT kit is used to determine acidity concentrations. Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 10-4,000 mg/L. Contains approximately 100 tests.

Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration

Rugged titration device is ideal for field use, with no glass burets required

Προδιαγραφές

Αριθμός δοκιμών :	100
Εύρος μέτρησης :	10 - 4,000 mg/L CaCO ₃
Κατηγορία:	Ψηφιακός τιτλοδότης
Μικρότερα βήματα ανάλυσης :	0.1-1.10
Μοντέλο σετ μέτρησης :	AC-DT
Όνομα μεθόδου :	Digital Titrator/Sodium hydroxide
Παράμετρος :	Acidity - (P or Total) as CaCO ₃
Τύπος θήκης :	F
Υποσημείωση :	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. 1EPA-accepted for wastewater. Check with your regulatory agency.