



Test kit, carbon dioxide, model CA-DT

Αρ. προϊόντος:

2064100

Επικίνδυνο

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής



Utilizing the Hach Digital Titrator, the Model CA-DT kit is used to determine carbon dioxide concentrations. Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 10.1-1,000 mg/L. Contains approximately 50 tests.

Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration

Rugged titration device is ideal for field use, with no glass burets required

Προδιαγραφές

Smallest increment:

0.1, 2.0

Αριθμός δοκιμών :

50

Εύρος μέτρησης :

10 - 1,000 mg/L CO₂

Κατηγορία:

Ψηφιακός τιτλοδότης

Μοντέλο σετ μέτρησης :

CA-DT

Όνομα μεθόδου :

Digital Titrator/

Sodium hydroxide

Phenolphthalein

Παράμετρος :

Carbon Dioxide - As CO₂

Τύπος θήκης :

F

Υποσημείωση :

*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. 1EPA-accepted for wastewater. Check with your regulatory agency.