



Bag, sterile, WhirlPak w/ dechlorinating agent, 100 mL, 25/pk

Αρ. προϊόντος:

2075325

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής

Specially designed for chlorinated water samples (Complies with USEPA and Standard Methods criteria). Sterile (certified), disposable, transparent polyethylene bags with trusted "Whirl-Pak" closure. Contains non-nutritive pill with 10 mg sodium thiosulfate (will dechlorinate up to 10 mg/ Cl_2 in 100 mL sample) Four ounce and 100 mL fill-lines. "Write-on" strip accepts waterproof ink.

Προδιαγραφές

Όγκος/Μέγεθος πακέτου :

25

Χωρητικότητα :

177 mL