



Be Right™



Rack, drying

Αρ. προϊόντος:

2089000

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής

Epoxy-coated, steel-wire drying rack for test tubes, vials, flasks, beakers, and other glassware. Three rows of hair-pin looped supports to hold large-diameter necks. Fourteen loops total. Top row has 12 straight-pin for smaller tubes and glassware. Polystyrene tray fits into bottom rails to catch draining water. Dimensions: 43 x 17 x 48 cm (17 x 6.7 x 19"). Free-standing or wall mount option.

Προδιαγραφές

Βάθος :	48 cm
Μήκος :	430 mm
Όγκος/Μέγεθος πακέτου :	1 each