



Test kit, copper (free & total) color disc, model CU-6

Αρ. προϊόντος:

2194100

Επικίνδυνο

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής



Copper may exist in natural waters and effluents as a soluble salt or as suspended solids. A small amount is essential for plants and animals, and concentrations exceeding 0.1 mg/L are useful for controlling algae and plankton growth. Quantities ranging from 0.02 to 0.1 mg/L are toxic for some fish, so its use for treating fish disease requires careful monitoring. This kit is USEPA-accepted for waste water and potable water and approved by some states for potable water under specific conditions (Check with your regulatory agency). Range: 0.1-5 mg/L; Contains 50 Free and 50 Total tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Προδιαγραφές

Parameter:	Copper (free & total)
Smallest increment:	0.1
Αριθμός δοκιμών :	50 (Free), 50 (Total)
Εύρος μέτρησης :	0 - 4 mg/L Cu
Κατηγορία:	Χρωματικός δίσκος
Μοντέλο :	CU-6
Μοντέλο σετ μέτρησης :	CU-6
Όνομα μεθόδου :	Color-disc/Bicinchoninate, Hydrosulfite reduction
Παράμετρος :	Copper - As free & total Cu [#]
Τύπος θήκης :	B