



Filter, TSS glass fiber, pore size 1.5 μm , diameter 20 mm, 100/pk

Αρ. προϊόντος: 2198720

EUR Τιμή: Επικοινωνήστε μαζί μας
Δεν εμφανίζεται η ημερομηνία αποστολής

TSS Filters are recommended for use in Standard Methods 2540 C and 2540 D. Optimized surface characteristics of this binderless borosilicate glass microfiber filter eliminates pin-holing and fiber-shedding to ensure consistent, accurate results. Fast flow rate, fine porosity and 1.5 μm particle retention makes the TSS Filter the preferred choice of environmental chemists for total suspended solids analysis of water and wastewater.

Προδιαγραφές

Διάμετρος :	20 mm
Μέγεθος πόρου :	1.5 μm
Όγκος/Μέγεθος πακέτου :	100 /pk
Υλικό :	Glass Microfiber