



## Test kit, general boiler water, drop count titration, model AL-94 MG-L

Αρ. προϊόντος:

222403

Επικίνδυνο

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής



Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. Used for controlling boiler and cooling water treatment and checking corrosion, this test kit is an economical solution for the small boiler operator who needs to monitor alkalinity, chloride, hardness, iron, and pH.

**Uses simple, economical drop count titration and color disc methods**

**Kit contains all required reagents and apparatus to perform analysis in the field**

**Uses drop count titration and color disc methods**

---

### Προδιαγραφές

Εφαρμογή:	Boiler Water
Κατηγορία:	Χρωματικός δίσκος
Μοντέλο :	AL-94 MG-L
Παράμετρος :	Drop Count Titration