



SA filter, cone, diameter 110 mm, soil analysis grade, 500/pk

Αρ. προϊόντος:

2263356

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής

Designed for sandy soils, SA filters are conically folded and have been optimized for use with the common extraction solutions used in soil nutrient analysis (Mehlich 1, Mehlich 3, Morgan, Modified Morgan, Bray, and Olsen). SA Filters provide the best combination of retention rating, filtration speed and purity requirements for soil nutrient testing (P, K, Ca, Mg). Medium filter retention.

Προδιαγραφές

Διάμετρος :	110 mm
Όγκος/Μέγεθος πακέτου :	50 /pk
Υλικό :	Cellulose