



Test kit, oxygen scavenger, color disc, model OS-1, 50 tests

Αρ. προϊόντος:

2349500

Επικίνδυνο

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής



Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. Utilizing the color disc method, this kit contains 50 tests.

Simple and fast color disc method

Measures common oxygen scavengers: DEHA, hydroquinone, ISA, and MEKO

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Προδιαγραφές

Αριθμός δοκιμών :	50
Εύρος μέτρησης :	0 - 52 µg/L DEHA (and other Oxygen Scavengers)
Εύρος μέτρησης (2):	0 - 1360 µg/L DEHA (and other Oxygen Scavengers)
Κατηγορία:	Χρωματικός δίσκος
Μοντέλο :	OS-1
Μοντέλο σετ μέτρησης :	OS-1
Όνομα μεθόδου :	Color disc/Iron reduction
Παράμετρος :	Oxygen Scavenger - As DEHA, hydroquinone, erythorbic acid, and methylethylketoxime
Τύπος θήκης :	R