



Be Right™



Test kit, boiler feed & scale, model BSC-1, drop count titration

Αρ. προϊόντος:

2350300

Επικίνδυνο

EUR Τιμή:

Επικοινωνήστε μαζί μας

Τηλεφωνήστε για την ημερομηνία αποστολής



Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. An economical test kit for monitoring scale formation in small industrial boilers, the BSC-1 Boiler Feed and Scale Test Kit uses drop count titration methods to determine levels of alkalinity, chloride, and hardness.

Προδιαγραφές

Εφαρμογή:

Boiler Water

Κατηγορία:

Ογκομετρική τιτλοδότηση