



Test kit, nitrogen, ammonia, model NI-SA

Αρ. προϊόντος:

2428700

Επικίνδυνο

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής



Designed for testing ammonia in freshwater or seawater, this kit's color disc method offers low-cost, fast, reliable results in an easy-to-carry kit.

Designed for testing ammonia in freshwater or seawater

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Προδιαγραφές

Smallest increment:	0.1
Αριθμός δοκιμών :	100
Εύρος μέτρησης :	0 - 2.4 mg/L NH ₃ -N
Κατηγορία:	Χρωματικός δίσκος
Μοντέλο σετ μέτρησης :	NI - SA
Όνομα μεθόδου :	Color disc/Salicylate
Παράμετρος :	Nitrogen, Ammonia, mid range - For freshwater or seawater**
Τύπος θήκης :	D