



## Sampler, water core, "sludge judge" large capacity, 1-1/4" diameter

Αρ. προϊόντος:

2506102

EUR Τιμή:  
Διαθέσιμο

Επικοινωνήστε μαζί μας

Accurate readings can be obtained on suspended solids up to a depth of 4.5 m (15'). Comes in three 1.5 m (5') sections with screw-type fittings. Designed for sampling solids in wastewater and food processing applications. Standard Sludge Judge features 3/4" diameter plastic tubing, while the large capacity Sludge Judge features 1-1/4" diameter plastic tubing for double the capacity.

---

### Προδιαγραφές

Διάμετρος :	31.75 mm
Όγκος/Μέγεθος πακέτου :	1 each