



Syringe filter, non-sterile PDVF, 33 mm, 0.45 um, 50/pk

Αρ. προϊόντος:

2513603

EUR Τιμή:

Επικοινωνήστε μαζί μας

Αποστολή εντός 2 εβδομάδων

Hach syringe filters offer an efficient way to remove particulates from samples prior to analysis. Filters feature 20% more surface area than standard 25 mm filters, which increases flow rates and capacity. PVDF (polyvinylidene fluoride) filters are ideal for solutions containing proteins. Standard Female Luer-Lok inlet with standard male slip Luer outlet. Hydrophilic. Diameter: 33mm. Pore Size: 0.45um. 50/PK.

Προδιαγραφές

Διάμετρος :	33 mm
Μέγεθος πόρου :	0.45 μm
Μονάδες μέτρησης:	50/pk
Όγκος/Μέγεθος πακέτου :	50 /pk
Υλικό :	PVDF