



## Parafilm M Laboratory Film

Αρ. προϊόντος:

251764

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής

Seals labware quickly, effectively. Film prevents sample evaporation and contamination; protects anhydrous materials. Can be folded repeatedly and twisted around sharp surfaces without rupturing. Stretched lengthwise at 21 °C (71 °F), Parafilm M increases approximately 200%. Becomes more flexible, softer, stronger as temperature increases. Resistant to alcohol, medium concentration acid, and alkalies. Width: 10 cm (4"); 38 m (125') roll.

---

### Προδιαγραφές

Όγκος/Μέγεθος πακέτου :

1 each