



Test Kit, AccuVac® Color Disc, Ozone, 0.05-1.50 mg/L

Αρ. προϊόντος:

2518050

Επικίνδυνο

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής



This ozone test kit uses the fast, accurate visual method for determining ozone concentrations. Use to control disinfectant levels in water, wastewater, and industrial water. Kit contains 25 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

AccuVac® vials mean no mess, no measuring

Indigo trisulfonate method masks interferences from chlorine and other oxidizing agents and makes ozone analysis easier

Rugged test kit contains all required reagents and apparatus to perform analysis in the field

Προδιαγραφές

Parameter:	Ozone
Smallest increment:	0.05
Αριθμός δοκιμών :	25
Εύρος μέτρησης :	0 - 1.50 mg/L O ₃
Κατηγορία:	Accuvac
Μοντέλο :	NA
Μοντέλο σετ μέτρησης :	AccuVac® Kit
Όνομα μεθόδου :	Color disc/ Indigo trisulfonate
Παράμετρος :	Ozone - As O ₃
Τυπικό εύρος μέτρησης:	0 - 1.50 mg/L O ₃
Τύπος θήκης :	G