



Bottle rack, plastic-coated, 15 compartments

Αρ. προϊόντος:

2551800

ΚΑΤΑΡΤΗΜΕΝΟ ΣΤΟΙΧΕΙΟ

Αυτό το στοιχείο δεν είναι πλέον διαθέσιμο.

Plastic-coated wire rack provides corrosion resistant transport system for sample bottles and bags. Racks feature 15 or 20 compartments with double-grid bottom to prevent loss of samples during transport. 15 compartment rack is ideal for Whirl-Pak® sample bags (4 oz). 20 compartment rack fits most bottles associated with standard drinking water sampling. Autoclavable.

Προδιαγραφές

Ειδικά χαρακτηριστικά :	Rack, 15 Bottle Compartments
Όγκος/Μέγεθος πακέτου :	1