



Test kit, surface water

Αρ. προϊόντος:

2559800

Επικίνδυνο

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής



The ideal tool for source water monitoring and wastewater treatment plants, this kit uses drop count titration and color disc methods to provide economical, easy-to-use tests for ammonia, chlorine, pH, nitrate, dissolved oxygen, phosphorus, and temperature. The rugged, lightweight kit includes the convenient pH Pocket Pal Tester and all necessary apparatus and reagents for testing.

Ideal tool for both industrial and municipal water treatment plants

Uses drop count titration and color disc methods to provide economical, easy-to-use tests for ammonia, chlorine, pH, titrate, dissolved oxygen, phosphorus, and temperature

Rugged, lightweight kit includes the convenient pH Pocket Pal[™] Tester and all necessary apparatus and reagents

Προδιαγραφές

Αριθμός δοκιμών :	100
Εφαρμογή:	Surface Water
Κατηγορία:	Χρωματικός δίσκος