



Be Right™



Sample transport kit

Αρ. προϊόντος:

2568700

EUR Τιμή:

Επικοινωνήστε μαζί μας

Τηλεφωνήστε για την ημερομηνία αποστολής

Contains everything needed for collecting field samples, keeping them cold during transport. Includes Rubbermaid Cooler, plastic-coated rack with 15 compartments for sampling bags, 100 Whirl-Pak sampling bags (177 mL) with dechlorinating agent, one leak-proof refrigerant pack. Dimensions: 42.4 x 28.4 x 58.6 cm (16.7 x 11.2 x 23.1")