



Cylinder, graduated, 1000 mL, w/ handles

Αρ. προϊόντος:

2612953

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής

One-piece polypropylene construction with a molded in handle top ensure a secure, solid grip and superior control for safer, more convenient pouring. The “Total Contained, Total Delivered” double scale is accurate and easy to read (cylinder will not create a meniscus), while the hexagonal base prevents tipping or rolling. Hand opening is 89x 29 mm (3 x 11/8"). Steam autoclavable upto 121°C (250°F).

Προδιαγραφές

Μέθοδος μέτρησης :	To Deliver/To Contain
Όγκος/Μέγεθος πακέτου :	1 each
Πιστοποιήσεις :	None
Σώμα υλικού:	Polypropylene
Τύπος :	Διπλό μετρικό
Χωρητικότητα :	1000 mL