



Bottle, wash, diamond alert, 500 mL, Variety 6/pk

Αρ. προϊόντος:

2613049

EUR Τιμή:

Επικοινωνήστε μαζί μας

Αποστολή εντός 2 εβδομάδων

Wide-mouthed and color-coded wash bottles provide safe and positive identification of common laboratory solutions. Standard Diamond Alert symbols indicate type of hazard contents may represent. Polyethylene bottles with color-coded polypropylene screw caps and spray-tip assembly. Variety 6-pack includes one bottle each for water, acetone, methanol, isopropanol, toluene, and ethanol.

Προδιαγραφές

Stem:	Blue Cap/Stem
Μονάδα:	6 pk
Όγκος/Μέγεθος πακέτου :	6 pk
Υλικό :	LDPE
Χρώμα :	Transparent
Χωρητικότητα :	500 mL