



## Sample vial, 120 mL w/ thiosulfate & custody tie, 200/pk

Αρ. προϊόντος:

2634001

EUR Τιμή:

Επικοινωνήστε μαζί μας

Αποστολή εντός 2 εβδομάδων

### Sample Vial with Thiosulfate and Custody Tie

120 mL sterile, hinged-cap vials meet all EPA requirements for microbiology sampling. Collection, incubation, and accurate reading can all be done in the vial without a transfer step. 100 mL raised fill line. Polypropylene. ( **Please Note:** These sample vials are not autoclaveable)

### Προδιαγραφές

|                         |                                 |
|-------------------------|---------------------------------|
| Διάμετρος :             | 43mm                            |
| Ειδικά χαρακτηριστικά : | Sample Vial, with Thiosulfate   |
| Μήκος :                 | 93mm                            |
| Μονάδα:                 | 200 /pk                         |
| Όγκος/Μέγεθος πακέτου : | 200 /pk                         |
| Χωρητικότητα :          | fill line at 100 mL (+/-2.5 mL) |