



Cylinder, graduated, mixing, certified, 25 mL

Αρ. προϊόντος:

2636340

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής

This certified, Class A glass mixing cylinder is calibrated “to contain” according to ASTM E542 & E1272. Serialized and traceable to NIST standards (includes Certificate of Identification and Capacity). Single metric scale in durable white markings with funnel-top. Glass penny head stopper included.

Προδιαγραφές

Μέθοδος μέτρησης :	To Contain
Όγκος/Μέγεθος πακέτου :	1 each
Πιστοποιήσεις :	Class A
Σώμα υλικού:	Glass
Τύπος :	Μονό μετρικό, Ανάμιξη
Χωρητικότητα :	25 mL