



## Cylinder, graduated, mixing, certified, 50 mL

Αρ. προϊόντος:

2636341

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής

This certified, Class A glass mixing cylinder is calibrated “to contain” according to ASTM E542 & E1272. Serialized and traceable to NIST standards (includes Certificate of Identification and Capacity). Single metric scale in durable white markings with funnel-top. Glass penny head stopper included.

---

### Προδιαγραφές

|                         |                       |
|-------------------------|-----------------------|
| Μέθοδος μέτρησης :      | To Contain            |
| Όγκος/Μέγεθος πακέτου : | 1 each                |
| Πιστοποιήσεις :         | Class A               |
| Σώμα υλικού:            | Glass                 |
| Τύπος :                 | Μονό μετρικό, Ανάμιξη |
| Χωρητικότητα :          | 50 mL                 |