



Filter Unit, Type A

Αρ. προϊόντος:

2656600

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής

Certified for water quality work, this sterile, disposable, polystyrene filter unit contains a 47 mm, 0.45 μ m cellulose nitrate membrane for superior recovery and growth of microorganisms. It offers good chemical resistance to dilute acids, aromatic and aliphatic hydrocarbons, and conforms to American Public Health Association (APHA) and Standard Methods and U.S. EPA water quality testing.

Upper receivers are graduated from 20 to 115 mL. After filtration, upper chamber can be easily twisted off of receiver.

Προδιαγραφές

Διάμετρος :

47 mm

Όγκος/Μέγεθος πακέτου :

12 per pack