



Vacuum filter microfunnel, pk/50.

Αρ. προϊόντος:

2831500

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής

Feature a completely assembled filtration unit: Funnel, filter, nutrient pad and petri dish. Filter: GN-6 Metrical with 0.45 µm pore size and 47 mm diameter; white, gridded surface. Funnel: Polypropylene Petri Dish: polystyrene lid and polypropylene base Nutrient Pad: Cellulose. Fits 2831500, 2846100 and 2846150 manifolds. Sterile and individually bagged.

Προδιαγραφές

Διάμετρος :	47 mm
Μέγεθος πόρου :	0.45 µm
Όγκος/Μέγεθος πακέτου :	50
Υλικό :	Mixed-Cellulose Ester