



Filter membrane, EZ-Pak, 0.45 micron, 47 mm, white, gridded, 600/pk (millipore)

Αρ. προϊόντος: 2916800

EUR Τιμή: Επικοινωνήστε μαζί μας
Δεν εμφανίζεται η ημερομηνία αποστολής

Millipore EZ Pack Membranes, 0.45 μ m pore size, 47mm

Dispenser provides superior convenience and assures sterility. Unique one press on the lever and a sterile membrane is presented (horizontally and inverted) for efficient and convenient assembly. The EZ Pak membranes come as a band of membranes that easily loads into the dispenser. Packaging materials are wound into the dispenser ensuring a clear operating area and quick change of bands. The membranes are Type HA, which signifies compliance with Standard Methods criteria.

Προδιαγραφές

Διάμετρος :	47 mm mm
Μέγεθος πόρου :	0.45 μ m
Όγκος/Μέγεθος πακέτου :	600
Υλικό :	Mixed-Cellulose Ester