



Bottle, storage, autoclavable, wide mouth, 125 mL pk/12

Αρ. προϊόντος:

2939103

EUR Τιμή:

Επικοινωνήστε μαζί μας

Δεν εμφανίζεται η ημερομηνία αποστολής

Polyethylene bottles suitable for sample collection to chemical storage. Matched polypropylene caps have a valve-seal design for a perfect, leakproof fit. A protruding ring under the cap holds shrink-wrapping securely in place. Autoclavable 121°C (Do not autoclave with caps in place). Available in Narrow Mouth or Wide Mouth.

Προδιαγραφές

Body Material:	Plastic
Package size UOM:	12pk
Όγκος/Μέγεθος πακέτου :	12 /pk
Υλικό :	PP
Χωρητικότητα :	125 mL