



Disposable BOD bottles, 300 mL, 100/cs

Αρ. προϊόντος:

2943100

EUR Τιμή:

Επικοινωνήστε μαζί μας

Αποστολή εντός 2 εβδομάδων

Disposable BOD Bottle

Case of 100 disposable BOD bottles **without stoppers**. Plastic bottle is EPA approved for BOD measurement.

Single bottles available under article number 2943101. Plastic stoppers (25 pack) available under article number 2943900. Glass stoppers (quantity of 1) available under article number 190901.

Fresh bottle, ready for use

No cleaning required

Unbreakable and recyclable

Oxygen will not diffuse in or out of sample

EPA compliant

Προδιαγραφές

Body Material:	Plastic
Αριθμός μεθόδου :	8043
Όγκος/Μέγεθος πακέτου :	100 100/pk
Υλικό :	PET
Χωρητικότητα :	300 mL