



Xenosep consumables kit for EPA method 1664A testing

Αρ. προϊόντος:

2943800

EUR Τιμή:

Επικοινωνήστε μαζί μας

Αποστολή εντός 2 εβδομάδων

Order the Xenosep Consumables Kit for a quick and easy way to replace your most used components.

Consumables kit includes: 24 Xenosep SPE Filters, 47mm diameter, 24 Transfer Pipets, 15 mL capacity per pipet, 24 Inline Columns with Sodium Sulfate, 10g per column.

Προδιαγραφές

Όγκος/Μέγεθος πακέτου :

1 each