



MultiTest kit case (142 x 312 x 185 mm), blue polypropylene

Αρ. προϊόντος:

4660400

EUR Τιμή:

Επικοινωνήστε μαζί μας

Διαθέσιμο

Durable blue polypropylene carrying case allows easy transport of items for field testing.

Recessed handle in lid provides easy stacking

Contents are carried in upright position

Includes recessed flat-panel areas for labeling

Προδιαγραφές

Διαστάσεις (ΥxΠxB): 142 mm x 312 mm x 185 mm

Τύπος θήκης : F

Ύψος : 4.5 in