



QL filters, pore size 8 - 12 μm , diameter 150 mm, circular, 100/pk, qualitative

Αρ. προϊόντος: 50658

EUR Τιμή: Επικοινωνήστε μαζί μας
Δεν εμφανίζεται η ημερομηνία αποστολής

QL Filters are recommended for use in various environmental and wet chemistry reference methods. Extremely pure cellulose filters of nearly 100% alpha content are specially processed to yield negligible extractable levels in the determination and identification of unknown sample materials.

Προδιαγραφές

Διάμετρος :	150 mm
Όγκος/Μέγεθος πακέτου :	100 /pk
Υλικό :	Cellulose