



## QL filters, pore size 8 - 12 $\mu\text{m}$ , diameter 185 mm, circular, 100/pk, qualitative

Αρ. προϊόντος: 50659

EUR Τιμή: Επικοινωνήστε μαζί μας  
Δεν εμφανίζεται η ημερομηνία αποστολής

QL Filters are recommended for use in various environmental and wet chemistry reference methods. Extremely pure cellulose filters of nearly 100% alpha content are specially processed to yield negligible extractable levels in the determination and identification of unknown sample materials.

---

### Προδιαγραφές

|                         |           |
|-------------------------|-----------|
| Διάμετρος :             | 185 mm    |
| Όγκος/Μέγεθος πακέτου : | 100 /pk   |
| Υλικό :                 | Cellulose |