



Be Right™



Τεστ σε φιαλίδια LATON

Αρ. προϊόντος:

APC138

EUR Τιμή:

Επικοινωνήστε μαζί μας

Διαθέσιμο

Επικίνδυνο



Hach APC cuvette tests: Minimum amounts of reagents, maximum safety.

100 cuvettes for 50 Total N evaluations designed for use with the AP3900 Laboratory Robot.

Measuring range: 1 - 16 mg/L TN_b.

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in 4 packs with 25 cuvettes each. Officially approved for consent limits. Environmentally friendly due to less chemicals.

Precise and reliable measurement values

Maximum safety for users, thanks to the closed system, low amounts of reagents and included recycling/disposal service

Convenient and error-free dosing of the reagents

Barcode label for automatic test recognition in the photometer, including batch number, documented shelf life, Hach's Truecal functionality and COA via RFID technology

Differentiation between tests and measuring ranges by means of colour-coding

Προδιαγραφές

Αριθμός δοκιμών :	50
Εύρος μέτρησης :	1 - 16 mg/L TN _b
Ημερομηνία Λήξης :	18 μήνες από την ημερομηνία παραγωγής
Κατηγορία:	APC
Όνομα μεθόδου :	EN ISO 11905-1
Παράμετρος :	Ολικό άζωτο (Laton)
Συνθήκες Αποθήκευσης:	15 - 25 °C