



Be Right™



USB + Power Module for 2100Q/2100Q IS & DR1900

Αρ. προϊόντος:

LZV813.99.00001

EUR Τιμή:

Επικοινωνήστε μαζί μας

Διαθέσιμο

The USB + Power Module is an optional accessory for the 2100Q & 2100Q IS portable turbidimeters, and the portable spectrophotometer DR1900. This module adds further functionality to the instrument allowing the user to:

Transfer data to a PC in XML format

Transfer data to the Citizen PD-24 printer (2960100)

Operate the instrument on line power (100-240Vac, 50/60 HZ) and Charge NiMH batteries (2971304)

Perform firmware updates

Quickly enter Sample and User IDs with the optional barcode scanner (LZV566)

Προδιαγραφές

Output Voltage:

9V

Είσοδος ρεύματος :

100 - 240 VAC, 50/60 Hz

Σκοπός προμήθειας :

Includes module, module cover, power supply with four plug adapters, standard USB cable with Mini-USB connector, special USB cable adapter, and manual.