



## Titralab AS1000 Series Sample Changer, 20 positions, 90 mL beakers, for one probe

Αρ. προϊόντος: AS1000.99.20090

EUR Τιμή: Επικοινωνήστε μαζί μας

Αποστολή εντός 1 εβδομάδας

### The AS1000 Sample Changer is the perfect laboratory assistant for automating analyses on sample series.

The Titralab AS1000 series is an automatic sample changer used in analytical laboratories in connection with AT1000 Series titrators. The AS1000 is available in four configurations of different tray capacity that holds the sample beakers. The rotating module turns the tray to change the sample position while the lifting module moves sensors and tubes up and down into each sample. The AS1000 ensures automatic stirring of samples thanks to its integrated magnetic stirrer, and rinsing of probes between each sample analysed thanks to dedicated rinse positions. The AT1000 Series/AS1000 Series setup is entirely PC-controlled with the TM1000 software to automate sample analysis series. It archives results, covering a wide variety of aqueous sample applications such as pH/Conductivity/Alkalinity and Hardness in water, Chlorides content in food products, pH/Acidity in wine and juices, Free & Total SO<sub>2</sub> in wine.

#### Run sample analyses in total autonomy, without user intervention

The rotating module turns the tray to change the sample position while the lifting module moves sensors and tubes up and down into each sample. Single-radius sample trays for up to 30 samples.

#### Complete automation to free up your time

AS1000 Sample Changer multi-parameter automation solutions release operators' time from cumbersome repetitive analyses for more valuable tasks that cannot be automated.

---

### Προδιαγραφές

Number of sample positions: 20

Εγγύηση : 24 μήνες

Τι περιέχει το καλάθι: AS1000 Sample Changer contains: Lifting module, Rotating module with carousel tray, support base with magnetic stirrer, probe & tubes holder, RS-232 cable with USB/Serial adapter to connect AS1000 to PC, Power supply, set of magnetic stir bars, set of PP beakers & TM1000 PC-Software with ethernet cables to connect AT1000 to local network or PC.

---

### Τι περιέχει το καλάθι

AS1000 Sample Changer contains: Lifting module, Rotating module with carousel tray, support base with magnetic stirrer, probe & tubes holder, RS-232 cable with USB/Serial adapter to connect AS1000 to PC, Power supply, set of magnetic stir bars, set of PP beakers & TM1000 PC-Software with ethernet cables to connect AT1000 to local network or PC.